CLASSIC TRINOCULAR



Objectives (Anti fungus coated)

- IOS Plan Achromatic 4X
- IOS Plan Achromatic 10X
- IOS Plan Achromatic 20X(S.L.)
- IOS Plan Achromatic 40X Oil(S.L.)
- IOS Plan Achromatic 100X Oil(S.L.)

Mechanical Body and Pre-Focusing

- Co-axial focusing system. PRE-FOCUSING LEVER & TENSION ADJUSTMENT RING.
- Reverse angle quintuple nosepiece revolves on ball bearing; comes with a Slot for polarizer.
- Double plate mechanical stage size 150mm X 140mm having low positioned co-axial controls.

Sub Stage Condenser

- Abbe condenser N.A.1.25 swing out filter holder and iris diaphragm.

Illumination System

- Brightness control knob: Kohler illumination system with 3W LED.

Standard Accessories

- Operating manual, dust cover, guarantee card and Styrofoam molded pack.

Trinocular Observation Head

- Siedentopf Observation head inclined at 30° and rotatable at 360°.
- Vertical focusable phototube for Digital Camera.
- Diopter adjustment ring on ocular tube (±5).
- Inter-pupillary distance from 48mm to 75mm.

Eyepiece (Anti fungus coated)

- Extra wide field 10X (paired) F.O.V. 22mm