CROWN TRINOCULAR



Objectives (Anti fungus coated)

- IOS Plan Achromatic 4X
- IOS Plan Achromatic 10X
- IOS Plan Achromatic 40X(S.L.)
- IOS Plan Achromatic 100X Oil(S.L.)

Mechanical Body and Pre-Focusing

- Co-axial focusing system. PRE-FOCUSING LEVER & TENSION ADJUSTMENT RING.
- Reverse angle quadruple nosepiece revolves on ball bearing.
- Double plate mechanical stage size 140mm X 132mm having low positioned co-axial controls On ball bearing guide ways.

Sub Stage Condenser

- Abbe condenser N.A.1.25 swing out filter holder, iris diaphragm.

Illumination System

- Illuminator systems with ultra high brightness 3W LED and Brightness control knob

Standard Accessories

- Operating manual, dust cover, guarantee card and Styrofoam molded pack.

Trinocular Observation Head

- Siedentopf Observation head inclined at 30° rotatable at 360°.
- Vertical focusable phototube for Digital Camera.
- Diopter adjustment ring on ocular tube (±5).
- Inter-pupillary distance from 48mm to 75mm.

Eyepiece (Anti fungus coated)

- Wide Field 10X (paired) F.O.V 22mm.