VICTORY TRINOCULAR



Mechanical Body and Pre-Focusing

- Quadruple nosepiece.
- PRE-FOCUSSING LEVER & TENSION ADJUSTMENT RING are on focusing mechanism.
- Hard coated mechanical stage is 180mmX165mm.Movement 50mmX40mm.

Illuminator

- 12V/20W Halogen lamp with spherical collector lens, continuous regulation of the illumination .

Condenser

- Long working distance 45mm.Condenser with N.A. 0.35. Field diaphragm, Pull out filter socket.

Filter

- Blue, Green, Yellow Glass Filter.

Standard Accessories

- Electric Cable, Fuse, Vinyle dust cover, Operating Manual, 12V/20W Halogen bulb, STYRO FOAM molded pack .

Trinocular Observation Head

- Siedentopf Observation head inclined at 45˚.
- Diopter adjustment ring on ocular tube.
- Inter-pupillary distance from 48mm to 75mm.

Eyepiece

Eyepieces EWF 10X (Paired) FOV 18mm.

Objectives (Anti fungus coated)

- LWD Plan 4X
- LWD Plan phase 10X
- LWD Plan phase 20X PH
- LWD Plan 40X