VICTORY PLUS (UPGRADED)



Objectives (Anti fungus coated)

- LWD Plan objective 4X
- LWD Plan objective phase 10X
- LWD Plan objective phase 20X
- LWD Plan objective phase 40X

Mechanical Body

- Quintuple nosepiece. PRE-FOCUSING LEVER & TENSION ADJUSTMENT RING are on focusing mechanism.
- Hard coated fixed mechanical stage is 250mm X 160mm, slider clip, extension plate and co-axial Movable stage can be assembled.

Illumination System

- 6v/30w Halogen Bulb. External transformer, intensity adjustable.

Sub Stage Condenser

- N.A-03 Koehler illumination condenser, with Phase-adjustable slide plate, long working 72 mm.

Filter

- 45mm green interference filters, color temperature change filter, infrared eliminate filter.

Standard Accessories

- Electric Cable, Fuse, Vinyle dust cover, Operating Manual, 12V/20W Halogen bulb, Styrofoam molded pack.

Trinocular Observation Head

- Siedentopf Observation head inclined at 45˚.
- Vertical phototube for micro photography.
- Diopter adjustment ring on ocular tube.
- Inter-pupillary distance from 54mm to 75mm.

Eyepiece

- PWF 10X (paired) F.O.V. 22mm and centering telescope.